

Ultrasonic Heating Cleaner|
Ultrasonic Heating Cleaning Machine
(BDT Series)



Product Description

BDT series ultrasonic cleaner has not only the cleaning function, but also the functions of extracting, emulsifying, accelerating dissolution, crushing and dispersing. Widely used in: machinery, electronics, plastics, instrumentation, environmental protection, medicine, packaging, military, aerospace, shipping, automotive and other industries manufacturing and maintenance cleaning; Cleaning of straws, suckers and utensils, degassing before chromatography, cleaning of medical instruments, medical materials and utensils; Cleaning of jewelry, jewelry, watches, precious metals, precious stones, COINS, glasses, etc.

Product Features

1. Large TFT screen Display, digital display.
2. Working time countdown display.
3. Room temperature - 80 °C temperature setting range
4. Real time display the actual temperature in the cleaning tank
5. 1-99 minutes total working time setting
6. Equipped with a special stainless steel basket, flat cover.
7. Made of high quality stainless steel
8. The operation procedure of the instrument adopts the single-chip microcomputer software

Technical Specification

Model	Tank Size (L*W*H)	Volume (L)	Frequency (KHZ)	Power (W)	Time (min)	Rack	Drainage
BDT-100DT	300·150·100	4.5	40	100	1-99	Yes	/
BDT-120DT	240·140·100	3	40	120	1-99	Yes	/
BDT-120DT	240·140·150	5	40	120	1-99	Yes	/
BDT-3200DT	300·155·150	6	40	180	1-99	Yes	Yes
BDT-5200DT	300·240·150	10	40	200	1-99	Yes	Yes
BDT-5200DT	300·240·150	10	40	250	1-99	Yes	Yes
BDT-5200DT	300·240·150	10	40	300	1-99	Yes	Yes
BDT-4200DT	300·240·200	14.4	40	400	1-99	Yes	Yes
BDT25-12DT	500·300·150	22.5	40	500	1-99	Yes	Yes
BDT25-12DT	500·300·150	22.5	40	600	1-99	Yes	Yes
BDT-1000DT	600·300·300	54	28	1000	1-99	/	Yes
BDT-1200DT	600·400·300	70	28	1200	1-99	Yes	Yes

BDT-1500DT	600-500-300	90	28	1500	1-99	/	Yes
BDT-1800DT	600-600-300	108	28	1800	1-99	/	Yes
BDT-2000DT	650-600-300	117	28	2000	1-99	/	Yes
BDT-3600DT	750-600-300	135	28	3600	1-99	/	Yes
BDT-6000DT	1200-650-350	273	28	6000	1-99	/	Yes
BDT-8000DT	1800-650-350	410	28	8000	1-99	/	Yes

Application

Widely used in: machinery, electronics, plastics, instrumentation, environmental protection, medicine, packaging, military, aerospace, shipping, automotive and other industries manufacturing and maintenance cleaning; Cleaning of straws, suckers and utensils, degassing before chromatography, cleaning of medical instruments, medical materials and utensils; Cleaning of jewelry, jewelry, watches, precious metals, precious stones, COINS, glasses, etc.