


Three Tanks Thermostatic Water Bath

(BESM Series| Ambient~300℃)



Product Description

BCXSC series three tanks thermostatic water bath is equipped with three groups of independent flume and three corresponding temperature control devices, which can set, control and display the temperature of each group of bath at different temperatures at the same time. It is specially used for precise constant temperature and auxiliary heating in universities, industrial and mining enterprises and scientific research units.

Product Features

1. Three sets of controllers and three sets of tanks can be used at the same time, which is convenient for reaction comparison and non-interference.
2. Wide temperature range: room temperature to 100 °C and room temperature to 200 °C.
3. PID temperature control technology, precision up to 0.3 °C
4. Compact and reasonable structure, save lab space

Application

Three tanks thermostatic water bath is specially used for precise constant temperature and auxiliary heating in universities, industrial and mining enterprises and scientific research units.

Technical Specification

Model	BESM-8D	BESM-8A	BESM-10D	BESM-10A
Temperature Range °C	ambient~100	ambient~100	ambient~100	ambient~100
Temperature fluctuation °C	±0.3	±0.3	±0.3	±0.3
Temperature Volatility °C	0.1°C	0.1°C	0.1°C	0.1°C
Volumn(L)	5L	5L	8L	8L
The Power Supply	220V 50HZ	220V 50HZ	220V 50HZ	220V 50HZ
Tank Size(mm) W*D*H	110×100×130 ×3 tanks	110×100×130 ×3 tanks	200×200×150× ×3 tanks	200×200×150× ×3 tanks
Shape Size(mm) W*D*H	480×245×320	480×245×320	800×350×350	800×350×350

Service

1. 12 months warranty for thermostatic water bath.
2. Technical support is also available timely.
3. Reply your enquiry within 24 working hours.
4. OEM is available if the quantity is large.

-
5. Package and shipping worldwide.